



Supports patient in an ideal flexed position during spinal and proctological operations. Thanks to its curved structure, it supports the chest and pelvis, relieves pressure in the abdominal area and thus ensures a high level of comfort for patient.

- Designed to meet all requirements with 60 degree radiolucent features.
- Posture control with improved surgical access for Laminectomy control for improved surgical access
- Easy to use crank system
- The design of the Radiolucent spine frame allows easy transport and storage

#### STANDARDS

**Classification** Class I

**EU Conformity** Declaration of conformity

#### PACKAGE DIMENSIONS (WxLxH)

620x870x400 mm - 19 kg

#### DIMENSIONS (WxLxH)

500x750x250 mm

## TECHNICAL DRAWING

