

DATA SHEET

Portable Suction Unit, 220V

Model : SU60.06 / SU60.06T



SU60.06



SU60.06T

Inspital portable suction units can be carried and operated easily due to practical design. Designed to be used in small operations and emergencies. Portable design suitable for rapid interventions.

STANDARDS

IEC 60601-1:2019

ISO 10079-1:2022

ISO 14971

MATERIAL SPECIFICATIONS

Body Material ABS plastic

PACKAGE DIMENSIONS

Main Unit 230x420x360 mm - 6 kg

Stand 320x320x750 mm - 3 kg

GENERAL FEATURES

Max. Vacuum (kPa/mmHg)	~83kPa / 620mmHg
Flow	25 l/min
Noise level	45 dBA
Pump type	Oil free
Power supply	220–240V, 50–60 Hz

AMBIENT PROPERTIES

Humidity	20% – 85%
Atmospheric Pressure	0.49 atm – 1 atm


STANDART ACCESSORIES

Collection Jar	1pc, 2lt, reusable, polycarbonate
Silicon suction hose	1pc, 2m
Hydrophobic filter	3pcs
Catheter jar	1pc
Ajutaj	1pc
Energy cable	1pc

Please Note

Vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

TECHNICAL DRAWING



REV.01/10.24