



- Made of AISI 304 stainless steel designed to be used by the surgery team for hand and nail cleaning before entering the sterile zone
- Photocell display screen that shows the temperature and duration of water alternately from each faucet
- Specially designed with round concave corners, easy to clean and antibacterial material considering the areas of use & user's requests and to be compatible with international hygienic norms
- Photocell and knee controlled water flow and knee controlled soap flow at each station
- Thermostatic mixing valve for cold and hot water

#### PACKAGE DIMENSIONS

<b>Dimensions</b>	800x920x1440 cm (single) - 68 kg
	800x1720x1440 cm (double) - 120 kg
	800x2300x1440 cm (triple) -160 kg

#### ELECTRICAL FEATURES

<b>Input Voltage</b>	24 V / 50-60 Hz
<b>Electrical Protection Class</b>	Class I

#### MATERIAL SPECIFICATIONS

<b>Body Material</b>	t:1mm stainless plate
<b>Sink Material</b>	t:1,2mm stainless plate
<b>Knee Contact Part</b>	Polyurethane

## INFRASTRUCTURE REQUIREMENTS

Hot water installation

Cold water installation

Grounded electrical socket

Water drainage

Model	L (mm)	W (mm)	H (mm)	Explanation
SK60.10	800	635	1180	Single, Water Photocell
SK60.20	1600	635	1180	Double, Water Photocell
SK60.30	2200	635	1180	Triple, Water Photocell
SK61.10	800	635	1180	Single, Water-Soap Photocell
SK61.20	1600	635	1180	Double, Water-Soap Photocell
SK61.30	2200	635	1180	Triple, Water-Soap Photocell



## TECHNICAL DRAWING

