

Instrument Table

- Made of AISI 304 stainless steel
- Designed for operating and intervention rooms
- Four pcs 75 mm castors two with brake and two electrically conductive
- Bumpers on both casters
- Curved corners

Mechanical properties

Yield 0.2%-Rp 0.2 :235
 Yield 1%-Rp 1 :265
 Strength tensile-Rm :550
 Elongation A50 :45

International Designations

EN 10088 Number :1.4301
 EN 10088 Name :X5CrNi18-10

Technical Drawing

Article Codes	L (mm)	W (mm)	H (mm)	Weight(kg)
AB00.21	900	650	850	15
AB00.23	800	600	850	14
AB00.22	700	550	850	13

