

Waste Bins and Kick Bucket

- Made of AISI 304 stainless steel
- Lid and carriage handle (AB90.01, AB90.02)
- Rubber coated metal bumper against impacts
- Tip resistant
- Swivel castors
- Optional MR compatibility(AB90.01)
- Four pcs 75 mm castors two with brake and two electrically conductive

Mechanical properties

Yield 0.2%-Rp 0.2 : 235
Yield 1%-Rp 1 : 265
Strength tensile-Rm : 550
Alongation A50 : 45

International Designations

EN 10088 Number : 1.4301
EN 10088 Name : X5CrNi18-10



Article Codes	L (mm)	H min (mm)	Explanation	Weight(kg)
AB90.01	350	350	20 L	3,2
AB90.02	530	530	50 L	4,2
AB90.03	400	200	No Lid	

Technical Drawing

