

Goepel leg holder, pair

OT60.08

General features:

Model	: OT60.08
Product explanation	: Geopel Leg Rest, pair
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Gross weight	: 7,6 kg
Packing dimensions	: 440x510x270mm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

Product Features:

Body material	: Ø18mm chrome plated steel rod
Supportmaterial	: Polyurethane



Technical Drawing

