

Lateral Support

General features:

Model	: OT60.10
Product explanation	: Lateral support, pair
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Grossweight	: 6 kg
Packing dimensions	: 32x45x22cm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)
Body material	: Ø12mm stainless stick
Support material	: Polyurethane

Product Features:

- Pair
- Including fixing clamps
- Extendible
- Rotatable and tiltable
- Adjustable height
- Made of stainless steel

Technical Drawing

