

**Rinsing Bowl, rectangular
OT60.11**

General features:

Model	: OT60.11
Product explanation	: Irrigation Bowl
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Bowlfeatures	: Stainless Bowl
Bowldimensions	: 325x265x90 mm
Grossweight	: 1.9 kg
Packing dimensions	: 350x400x140 mm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

Product Features

Bowlmaterial	: t:1mm stainless profile
--------------	---------------------------



TECHNICAL DRAWING

