

Rinsing Bowl, Rectangular with Drainage Tube OT60.12

Model : OT60.12
Grossweight : 5 KG
Packing Dimensions : 35x45x15cm

- 20 L capacity
- Compatible with irrigation bowls with drainage
- Tube hole on the cover
- 5 castors and rubber bumper
- Including carriage handle and cover
- Made of stainless steel



Rinsing Bowl, Round with Drainage Tube OT60.13

Model : OT60.13
Grossweight : 5 KG
Packing Dimensions : 35x45x15cm

- Including fixing clamp
- Including drainage tube
- 2 joints rotatable arms
- Adjustable height
- Ø260x80 mm
- Made of stainless steel



Waste Bin OT60.14

Model : OT60.14
Grossweight : 5 KG
Packing Dimensions : 60x40x60cm

- 20 L capacity
- Compatible with irrigation bowls with drainage
- Tube hole on the cover
- 5 castors and rubber bumper
- Including carriage handle and cover
- Made of stainless steel

