

Arm / Hand Surgical Table**General features:**

Model	: OT60.26
Product explanation	: Armsurgery plate
Material	: Stainless steel(AISI 304, DIN X5 CrNi 18-10)
Wheel features	: 2 pcs Ø72,5mm wheel
Grossweight	: 12 kg
Packingdimensions	: 61x76x26cm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

Product Features:

Body material	: 30x30mm stainless box profile
Matressmaterial	: Anti bacteria limitation coated sponge floor

**TECHNICAL DRAWING**