

Orthopedic Traction Set Trolley

Model : OT60.33

Grossweight : 10 Kg

Packing Dimensions : 60x90x96cm

- For the storage and easy installation of orthopedic traction set
- Large basket for smaller accessories
- Four castors, two with brake
- Made of stainless steel



Accessory Trolley

Model : OT60.32

Grossweight : 12 Kg

Packing Dimensions : 60x90x96cm

- For the storage of table accessories
- Four rails to hang all accessories
- Four castors, two with brake
- Made of stainless steel

