



Gel Prone Head Rest

General features:

Model	: OT60.57
Product Explanation	: Prone Head Support
Material	: Silicone gel
Gross Weight	: 5.5 Kg
Packing Dimensions	: 30x25x18cm
Product Documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

TECHNICAL DRAWING



Spinal Disc Pad

General features:

Model	: OT60.61
Product explanation	: Wilson Frame Set
Material	: 1 Piece of large cushion : 4 Pieces of gel pads
Gross Weight	: 7 Kg
Packing Dimensions	: 20x75x53cm
Product documentation	: CE
Classification	: Class I (MDD 93/42 EEC)

TECHNICAL DRAWING

