



Leg positioning boot facilitates the posture of the patients to meet your specific needs in bariatrics, gynecology, lithotomy, urology or long lasting operations.

- Safe and easily adjustable in various positions
- Adjustable clamps for stabilization, remains securely in position
- Special comfort boots for pressure relief of the calves thanks to soft padded insert
- Prevents injury in the area of the nervus peroneus

STANDARDS

Classification	Class I
EU Conformity	Declaration of conformity

PACKAGE DIMENSIONS (WxLxH)

1000x510x310 mm – 16 kg

TECHNICAL SPECIFICATIONS

Cantilever	Stainless steel (304)
Cantilever Length	100 cm
Applicable Patient Weight	Max. 180 kg
Vertical Adjustment Range of the Cantilever	+90 ° / -45 °

Leg boot pads are made of soft memory foam, which is breathable and prevents pressure sores during long-term surgery.

TECHNICAL DRAWING

