

Neurosurgery skull clamp & horseshoe head rest OT60.80

General features:

Description

- Neurosurgery head rest especially designed for use in neurosurgery operations.

| | |
|--------------------------|--|
| Model | : OT60.80 SET |
| Product explanation | : Neuro surgery head rest with pin |
| Material | : Stainless steel(AISI 304, DIN X5 CrNi 18-10) |
| Gross weight | : 19 kg |
| Packing dimensions | : 23cm x 62cm x 55cm |
| Product documentation | : CE |
| Classification | : Class I (MDD 93/42 EEC) |
| Body material | : Antistatic Coated Aluminum Casting |
| Pin Material | : Stainless steel |
| Horizontal movement rail | : Ø32x2mm Aluminium tube |



Product Features:

- Extremely durable and rigid construction from aluminum alloy
- 3 pcs fixing system with pins and completely adjustable skull clamps
- Adjustment mechanism from 4 different points
- Easily mountable horse shoe head rest accessory

A : X - 360 degree
: Y - 60 degree
: Z - 60 degree

Technical Drawing

(min-max)122 mm - 300 mm

