

IV Stand

Item Code:
SR10.70 Turquoise Color
SR10.71 Anthracite Color

Configuration

- Hook height can be adjusted between 1200-2100 mm
- 4 hooks to hang IV
- Infusion devices mountable body structure
- 2 options with plastic or chrome plated star base
- 5 pcs castors 2 with brake
- Ø 50 mm Castors

Article Codes	L min (mm)	H min (mm)	H max (mm)	Weight(kg)
SR10.70	590	1200	2100	5 kg
SR10.71	590	1200	2100	5 kg



Technical Drawing

