



- Thermostatic mixingvalve for hot and cold water
- Automatic flow control (water, disinfectant, soap) at each station, with all coming from the same faucet
- A drawer located at the front of the handwashing unit
- High surface quality, easy cleaning, and antibacterial properties of the sink thanks to Corian
- Adjustable flow times over time
- Splash-resistant sink design

MATERIAL SPECIFICATIONS

Body Material	t:1mm stainless plate
Sink Material	12 mm Corian
Knee Contact Part	12 mm Corian

ELECTRICAL FEATURES

Input Voltage 24 V / 50-60 Hz

Electrical Protection Class Class I



INFRASTRUCTURE REQUIREMENTS

Hot water installation

Cold water installation

Grounded electrical socket

Water drainage

Model	L (mm)	W (mm)	H (mm)	Explanation
SK61.10	720	625	1780	Single, Water-Soap Photocell
SK61.20	1420	640	1780	Double, Water-Soap Photocell

TECHNICAL DRAWING

