



- Made of AISI 304 stainless steel designed to be used by the surgery team for hand and nail cleaning before entering the sterile zone
- Photocell display screen that shows the temperature and duration of water alternately from each faucet
- Specially designed with round concave corners, easy to clean and antibacterial material considering the areas of use & user's requests and to be compatible with international hygienic norms
- Photocell and knee controlled water flow and knee controlled soap flow at each station
- Thermostatic mixing valve for cold and hot water

PACKAGE DIMENSIONS

Dimensions	800x920x1440 mm (single) - 68 kg
	800x1720x1440 mm (double) - 120 kg
	800x2300x1440 mm (triple) -160 kg

ELECTRICAL FEATURES

Input Voltage	24 V / 50-60 Hz
Electrical Protection Class	Class I

MATERIAL SPECIFICATIONS

Body Material	t:1mm stainless plate
Sink Material	t:1,2mm stainless plate
Knee Contact Part	Polyurethane

INFRASTRUCTURE REQUIREMENTS

Hot water installation

Cold water installation

Grounded electrical socket

Water drainage

Model	L (mm)	W (mm)	H (mm)	Explanation
SK60.10	800	635	1180	Single, 1 Water - 1 Soap Photocell
SK60.20	1600	635	1180	Double, 2 Water - 2 Soap Photocell
SK60.30	2200	635	1180	Triple, 3 Water - 2 Soap Photocell



TECHNICAL DRAWING

