INSPITAL
"innovation for hospital"

DATA SHEET

Surgical Suction Unit with Extractor, 2x3L Reusable Colletion Jar
Model : SU60.05



Designed to be used in delivery operations via Malmstrom-Throe method by its integrated pedal system.
Device has optional delivery assistance equipment as birth cups and cannula sets. Various liter collection jar and liner options are available.

STANDARTS

IEC 60601-1:2019

ISO 10079-1:2022

ISO 14971

MATERIAL SPECIFICATIONS

Body Material ABS plastic

PACKAGE DIMENSIONS

Main Unit	510x510x960 mm - 6 kg
Stand	330x450x210 mm - 2.3 kg

INSPITAL MEDICAL TECHNOLOGY GmbH

P: +49 2131 291 75 00
W: www.inspital.com

Mosel Str. 23, 41464 Neuss
Germany

INSPITAL
"innovation for hospital"

DATA SHEET

Surgical Suction Unit with Extractor, 2x3L Reusable Colletion Jar
Model : SU60.05

GENERAL FEATURES

Max. Vacuum (kPA/mmHg)	~90kPa / 620mmHg
Flow	60 l/min
Noise level	55 dBA
Pump type	Oil free
Power supply	220–240V, 50–60 Hz

AMBIENT PROPERTIES

Humidity	20% – 85%
Atmospheric Pressure	0.49 atm – 1 atm

STANDARD ACCESSORIES

Collection Jar	2 pc, 3lt, reusable, polycarbonate
Silicon suction hose	2 pc, 2m
Yankauer hose	2
Hydrophobic filter	7 pcs
Catheter jar	1 pc
Ajutaj	2 pc
Energy cable	1pc


OPTIONAL ACCESSORIES

Safety jar	

Please Note

Vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

TECHNICAL DRAWING



INSPITAL MEDICAL TECHNOLOGY GmbH

P: +49 2131 291 75 00
W: www.inspital.com

Mosel Str. 23, 41464 Neuss
Germany